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- 0,0 will appear flashing on the screen.
- Set the difference value using the + or – buttons. With each pressing of the buttons, the value will increase/decrease by 0,1°C.
- Confirm the settings with SET.
- To turn off the thermometer, press SET 5 times or wait 5 seconds until the thermometer turns off automatically.

14. CLEANING AND SAFETY PRECAUTIONS

- The sensor of the device is fragile. To ensure most precise measurements, it must be clean and undamaged.
- Clean the sensor with cotton swabs lightly moistened with alcohol. Point the sensor downwards to prevent it from damping.
- Clean the casing with a soft cloth moistened with water or alcohol.
- Do not submerge the thermometer in water or other liquid.
- Remove the batteries when the thermometer is not in use for an extended period of time.
- Do not submerge the device in water or expose it to direct sunlight.
- Do not throw or hit the device; doing so will result in damage.
- Protect the product from moisture, dust and other dirt.
- The thermometer was designed for practical use but it does not substitute for a medical examination.
- Use and store the thermometer in the conditions of operation (see technical specification).

15. TROUBLESHOOTING

Problem and reasons	Advice
Low display contrast.	Check/replace the batteries.
Incorrect temperature or large disparity between measurements.	Check/replace the batteries or check the distance between the measured object and the thermometer. During measurement, the thermometer should be located 5-15 cm away from the scanned surface.
Measurement result shows “Lo”, temperature is beyond the thermometer’s range.	Such situation may occur in the Body mode when the temperature is below 32°C, or in the Surface mode when the surface temperature is below 0°C.
Measurement result shows “HI”, temperature is beyond the thermometer’s range.	Such situation may occur in the Body mode when the temperature is above 42,9°C, or in the Surface mode when the surface temperature is above 100°C.
The thermometer does not work, it does not perform any measurements.	Check/replace the batteries; or the ambient temperature is outside the scope of 5°C – 40°C. Take the thermometer to a room with appropriate temperature conditions.

This symbol, placed on electric and electronic devices, means that the purchased product should not be disposed of along with ordinary household waste. It should be taken to an appropriate collection and recycling point of electric and electronic devices. Proper disposal and recycling of equipment helps prevent the potentially negative consequences for the environment and human health.

Measurement:	Normal temperature (°C)	Normal temperature (°F)
Rectal	36.6 to 38	97.8 to 100.4
Oral	35.5 to 37.5	95.9 to 99.5
Ear	35.8 to 38	96.4 to 100.4
Forehead	35.8 to 37.8	96.4 to 100.0

Human body temperatures change during the day, as well as under the influence of many external factors such as age, sex, skin type and thickness, etc.

11. CLOCK SETTINGS

The device can display the current room temperature and the current time. To enable this function:

- Turn on the thermometer using the ON/OFF button (trigger).
- Press and hold MODE for about 3 seconds; "TIME" will appear on the screen and the hour number will begin flashing.
- Set the hour using UP or DOWN,
- Press SET, the minute number will begin flashing.
- Set the minutes using UP or DOWN.
- Press SET to confirm the settings.

Every 10 seconds the display will change, showing the time (TIME) and the room temperature (ROOM) alternately.

To disable this function, reset the thermometer by removing the batteries.

12. MEMORY OF THE THERMOMETER

Measurement results are recorded automatically. The thermometer records the last 50 measurement results. When the memory is full, the oldest data will be replaced by the new ones.

How to access the memory of the thermometer:

- Turn on the thermometer using the ON/OFF button (trigger).
- Press UP; the screen will display "MEMO" and "1", the number of the latest measurement.
- Pressing UP further, you can check the recorded measurement results in sequence.

If the recorded measurement result applies to body temperature, inscription "Body" will be additionally displayed next to the inscription MEMO. If the recorded measurement result applies to surface temperature, the display will show the inscription "Surface temp".

- To clear the memory of the thermometer, press and hold UP. The measurement results will be replaced by "----".

13. THERMOMETER CALIBRATION

Calibration should be performed at least once per month and at every change of environment in which the measurement is performed. To calibrate the device, measure a patient's armpit temperature in the environment (room) where the measurements are to be performed, with a traditional non-mercury thermometer and with the 050400 thermometer. If the results differ, follow the advice below.

- Calculate the difference in readings between the traditional (non-mercury) thermometer and 050400.
- To correct the difference, turn on the thermometer using the main button, all segments of the display will appear on the screen; when you hear a sound, press SET; the LCD screen will display F1.

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- Aim the thermometer at the forehead or ear of the examined person (see the picture below for positioning), 5-15 cm from the body.

- Press the measurement trigger.
- When the measurement is finished, you will hear a sound. The measurement result will be displayed in 0,5 second.
- The thermometer will turn off automatically after 5 seconds.

d) Surface temperature measurement

- Turn on the thermometer with the ON/OFF button (trigger).
- Press MODE to switch between the body mode (Body) and the surface mode (Surface temp). Set to “Surface temp”.
- Bring the thermometer 5-15 cm close to the surface of a liquid etc. Do not put the thermometer into the liquid.
- Press the measurement trigger.
- When the measurement is finished, you will hear a sound. The measurement result will be displayed in 0,5 second.
- The thermometer will turn off automatically after 5 seconds.

9. LIGHT AND SOUND SIGNALIZATION OF HIGH TEMPERATURE

The light signalization applies only to the Body mode, i.e. body temperature measurement. In surface measurements (in the Surface mode), every measurement will be backlit in green.

If the fever alert is set at 38°C (default settings):


- green backlight – when the temperature taken is below 32°C - 37,5°C
- orange backlight and 3 short sounds – when the temperature taken shows values between 37,6°C and 38,0°C
- red backlight and 3 short sounds – when the temperature taken is between 38°C and 42,9°C .

Remarks:

- The fever alert may be set in the body temperature range between 32°C and 42,9°C.
- If the alert temperature is set above 38°C, backlight will turn red when the measurement result exceeds this value.
- If the fever alert temperature is set to 37,6°C or below, the device will not display the orange backlight.

10. GENERAL KNOWLEDGE ON THE BODY TEMPERATURE

Human body is an extremely complex system, and body temperature is a very important indicator of health conditions of the organism. You can measure the temperature of your forehead, ears or mouth on a supplementary basis, to check your health conditions. For different body parts, the normal value of temperature is different. The table below shows these relations.

- Close the lid.
- The thermometer is ready to begin a measurement.

7. SETTINGS MADE USING THE “SET” BUTTON

Using the SET button, you can programme the following functions:

- Thermometer calibration **F1**
- Temperature display unit selection **F2**
- Fever alert settings **F3**
- Sound on/off **F4**

Function programming:

- Turn on the thermometer with the ON/OFF button (trigger).
- Press SET, the screen will display F1 (described in point 13).
- Press SET again, the screen will display F2, then °C will appear flashing.
- Set the temperature unit – Celsius (°C) or Fahrenheit (°F) degrees – using the UP or DOWN button.
- Press SET again, the screen will display F3, and the alert temperature, set by default to 38 °C, will flash.
- Set the alert temperature using the UP or DOWN button. Above the programmed temperature, the backlight will change to red.
- Press SET again, the screen will display F4.
- Using the UP or DOWN button, enable (ON) or disable (OFF) the sound signal of the thermometer.
- Press SET to confirm the settings and turn off the thermometer.

8. MEASUREMENT

Remarks concerning measurement of the forehead temperature:

- Measurement should only be performed on an uncovered part of the forehead; any covering (e.g. with hair or hat) may result in wrong measurement.
- In case of temperature measurement in babies, who are usually wrapped in a blanket or clothing, uncover the neck for at least 1 minute.
- Do not measure the temperature immediately after removing sweat from the forehead.
- Bear in mind that the forehead surface temperature is different in different areas due to differences in blood flow intensity. Therefore, temperatures of the left and right temple may differ.

Remarks concerning measurement of the ear temperature:

- In order to obtain a precise measurement, clear the ear of the excess of ear wax.
- Temperature measured in the right ear may differ from the measurement in the left ear. To obtain a precise measurement result, always measure the temperature in the same ear.
- Do not use the thermometer in case of ear diseases or when the ear canal is moist.
- In case of hair in the ear, the displayed temperature may be lower than the actual.

c) Body temperature measurement

- Turn on the thermometer with the ON/OFF button (trigger).
- Press MODE to switch between the body temperature measurement mode (Body) and the surface scanning mode (Surface temp). Set to Body.

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Operation Manual
Professional touchless body thermometer
Electronic, infrared
050400

This operation manual constitutes a part of the product and should be stored in a manner allowing to acquaint oneself with its content in the future. It contains important information on the settings and operation of the device.

1. THE SET INCLUDES:

- An electronic thermometer with an infrared sensor
- Thermometer storage pouch
- Operation manual

2. FUNCTIONS

- Infrared body thermometer – Temperature is taken in 0,5 second.
- Touchless – measurement from 5-15 cm
- Memory: 50 measurements
- Measurement temperature in °C or °F
- Scanning function – to measure the surface temperature of baby milk, bathing water etc.
- Large and clear LCD screen
- 3 backlight colours
- Calibration possible

3. BEFORE USE:

- The product is intended only for the use described in the present manual.
- A measurement with this product cannot substitute for medical examination.
- Do not scan scars or damaged tissues.
- A precise measurement requires cleaning of the scanned surface (e.g. forehead) of make-up and other substances which could distort the measurement result.
- Do not use near sources of heat (e.g. fireplaces, ovens).
- The device should be used in an ambient temperature between 5°C and 40°C (0°F and 140°F). Do not use the thermometer outside this range.
- To obtain most precise measurement results, the thermometer and the examined person should be in the same room for about 15-20 minutes to balance their temperatures before the measurement.
- Maintenance activities which may be performed unaided include only battery replacement and cleaning.

4. PRINCIPLE OF OPERATION OF THE THERMOMETER

Every object with a temperature higher than absolute zero gives off a quantity of infrared energy, proportional to its temperature. An infrared thermometer treats a surface of a scanned object as opaque. The optical system of the thermometer directs and focuses the infrared radiation energy in the detector. Subsequently, the electronic circuits of the device process the signal from the detector to a reading shown on the screen.


The 050400 touchless thermometer is specially designed to measure body temperature (Body Mode) and surface temperature (Surface Mode), it can also scan temperature of rooms, food, bathing water, baby milk, etc. It also has a clock function. Due to the use of advanced infrared technologies, the device measures the temperature and analyzes the data in 0,5 second.

Caution: Depending on the skin type and thickness, the measurement result may be different for individual persons.

5. DEVICE DESCRIPTION


Buttons:

1. Infrared sensor
2. LCD screen
3. MODE button
4. DOWN button
5. UP button
6. SET button
7. ON/OFF button - trigger
8. Battery compartment

Display description:

1. Surface Mode – surface temperature measurement
2. Body Mode – body temperature measurement
3. Room – Room temperature
4. Time – clock
5. Measurement result
6. Temperature measurement icon – happy or sad
7. Battery low volume indicator
8. Speaker icon
9. Memo – memory
10. °C
11. °F
12. Number of a recorded measurement

6. INSTRUCTION OF USE

● **Battery installation**

- Open the battery compartment lid in the lower part of the thermometer.

- Insert 2 AAA batteries into the compartment.

Caution: To avoid problems of operation, pay special attention to the polarization of alkaline batteries during installation (installation of batteries with wrong polarization may cause permanent damage to the device).